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
RPS 300

Integrated lift with CAN BUS managing system

Original instruction translation

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English text of this document is the translation of the original manufacturer instructions. If further assistance is required, please refer to the original document or forward a request to **Mobility Networks Ltd**.

This handbook cannot be reproduce, stored by a filing system or transmitted to third parties, in whole or in part, in whatsoever form or medium, whitout **Mobility Networks Ltd** prior written permission.

Mobility Networks Ltd reserve the right to carry out products changes and improvements in order to enhance their quality, at any time and without notice, even il such modifications will not he reparted in this pubblication.

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SAFETY RULES AND SYMBOLS

Some precaution are listed in this manual in order to identify potentially dangerous situation in advance and the relevant instructions to avoid them, identified by the following symbols:

 WARNING	This symbol identifies the presence of instructions which need to be read and followed carefully in order to avoid potentially dangerous situations.
---	--

 DANGER	This symbol identifies the presence of essential information needed to avoid potentially dangerous situations that could cause physical injuries and/or damages to the equipment.
--	---

To ensure smooth and safe operation, it is necessary to know the exact procedure for the installation and service of the wheelchair lift.

WARNING

Carefully follow the instructions for the installation and service of the wheelchair lift within this manual. If the instructions within this manual have not been fully understood or further information is required, please contact Mobility Networks Ltd immediately. These indications non-fulfilment can cause serious damages, to the equipment and to the vehicles as well as serious injuries to people.

DANGER

Should it be necessary to raise the vehicle being equipped, please check that the characteristics of the lifting device or of the mobile bridge are compatible with the volume and weight of the actual vehicle otherwise there is a risk of serious damage to equipment and vehicles as well as serious injuries to people.

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1. PREFACE

1.1. Classification and general considerations

The wheelchair lift "RPS 300", manufactured by **Mobility Networks Ltd**, which is the object of this instruction manual, is in compliance with the safety requirements of 2006/42/CE directive. **RPS 300** wheelchair lift, as per the aforementioned directive, must be considered a 'machine' and therefore a declaration of conformity has been released, according to annex II of 2006/42/CE directive.

All the necessary instructions for the correct installation of the wheelchair lift are detailed for the installer together with a list of instructions for servicing procedures and for extraordinary maintenance.

WARNING

All installation and service procedures must be carried out in compliance with current health and safety laws.

1.2. Duties of installer / service

Installer / servicer must:

- carry out the installation in compliance with the instructions detailed within this manual;
- carry out the installation following the instructions of the vehicle manufacturer;
- fill and sign the installation conformity declaration;
- carry out the verification check procedure for the first commissioning following the instructions detailed in this manual;
- register the verification check procedure for the first commissioning in the "Use and Maintenance Manual";
- give to the final customer all the documentation required by the rules and regulations in force;
- carry out the inspections in compliance with the instructions detailed in this manual;
- register the inspections in the "Use and Maintenance Manual".

1.3. Lift Packing

For transportation, lift and its accessories are packed in a box made of wood and cardboard.

To move the packed lift a fork lift truck or a hoist (with proper capability) is advised.

All packing materials are ecological and recyclable. The packing materials should be disposed off correctly.

The local council or refuse department will be able to offer advice regarding disposal requirements.

DANGER

Cardboard and protective plastic sheets used for packing purposes can cause suffocation.

Packing mass is 180kg, accessories included (according to the model).

WARNING

During the unpacking process be very careful to not damage the contents.

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1.4. Delivery composition

The following is the materials list (and their quantity) that contain the packaging:

- | | |
|---|-------|
| • Wheelchair Lift RPS 300 | 1 pcs |
| • Electro-hydraulic control | 1 pcs |
| • Wired remote control panel 4 Buttons | 1 pcs |
| • Support of wired remote control panel | 1 pcs |
| • Lug Ø8 sect.25 A5-M8 | 6 pcs |
| • 90° rubber cup | 1 pcs |
| • Enable-key (080997) | 1 pcs |
| • Two Poles cable (2x0.75 black) | 10 mt |
| • Led-light red blink | 1 pcs |
| • Fuse mega 100 cod.F17105 | 1 pcs |
| • Self fuse mega cod.F17113 | 1 pcs |
| • Black corrugate sheath | 15m |

1.5. Documentation

Final user must receive with the "RPS 300" wheelchair lift:

- the original copy of the "Manufacturer's declaration of conformity";
- the installation "Conformity declaration" both filled and signed by the installer;
- "Use and Maintenance" manual filled by the installer;
- Spare Parts Catalogue (on request).

Mobility Networks Ltd is at your complete disposal for further clarification and instructions.

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1.6. Warranty

Mobility Networks Ltd assures the wheelchair lift for 24 months from the date of delivery. The warranty covers for defects concerning material quality and product manufacturing.

The warranty does not cover consumables and defects or failures as a result of installation procedure incorrectly carried out or not in compliance with the manufacturer's instructions.

The warranty becomes void in case of impacts caused by accidents and/or tampering carried out by personnel not authorized by **Mobility Networks Ltd**.

Mobility Networks Ltd refuses all responsibility for damages caused by:

- improper use of the electro-hydraulic lift;
- lift maximum load capability (300kg);
- a failure in carrying out "Use and Maintenance" manual instructions;
- a failure in carrying out maintenance operation as detailed in the "Use and Maintenance" manual;
- interventions or modifications to the lift without **Mobility Networks Ltd** authorization;
- lift safety devices removed or tampering;
- wrong manoeuvres when moving the hydraulic lift;
- presence of objects or people in lift working area.

Manufacture's regular inspections dates non-fulfilment will automatically make void the warranty.

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2. GENERAL INDICATIONS

2.1. Technical Description

RPS 300 model wheelchair lift (**Fig.1**) has been exclusively designed for lifting and lowering disabled people and/or only one operator onto/from the vehicle floor to which it is fixed.

Electro-hydraulic lift **RPS 300** is mainly composed by:

- main steel frame (cassette) strictly anchored to the vehicle frame with special brackets;
- sliding carriage, with arms and lifting cylinders connected to it, that allows platform to exit from the structure;
- loading platform covered with antiskid aluminium sheet;
- electro-hydraulic power unit, with electro-pump, oil tank, emergency and electric valves and hand pump.


A side installation are normally preferred.

On an horizontal level, platform is moved by an electric gear-motor with rack gear system; on a vertical level, it is moved by hydraulic cylinders with parallelogram kinematics.

Lift is fed by vehicle batteries by means of a spanner switch and it is managed by a electronic – unit.

The manual remote control panel (keypad) is provided with push-buttons allowing platform enter/exit and lifting/lowering.

In emergency conditions, platform can be controlled not only by remote control, but also by lever switches fixed on the electric control panel. In addition, even if electric installation failure has occurred, the lift can be used in safety conditions by means of pump and valve manuals and unlocking handle.

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2.2. Technical Data

The following are the main technical specifications of the wheelchair “**RPS 300**”:

- loading capacity: 300 Kg
- weight: 160 Kg
- electrical system: 12 Volts
- maximum lifting height: from cassette a load level, 400 mm about;
- maximum lowering height from cassette a floor: 320 mm about;
- electro-hydraulic power unit:
P = 14 Mpa; cc/giro = 0,45
- equivalent continuous acoustic pressure (A):
<70 dB (A);


Fig.2

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2.3. Safety Devices

The wheelchairs lift “RPS 300” they are equipped with devices of protection and safety:

- **Hydraulic installation general maximum pressure valve:**
operating in order to reduce the hydraulic oil pressure to values not higher than the ones of the working pressure (required for lifting the maximum load allowed).
- **25 A fuse;**
in order to protect the electronic unit.
- **100 A fuse on the battery:**
place on the positive cable of the battery and before the enable-key for protection entire power line.
- **External flap to stop and connect platform to the ground and Internal flap to stop and connect platform to vehicle loading deck:**
to make wheelchair going up/down the platform easier; they are provided with torsion springs in order to keep them vertical when platform is far from ground or vehicle loading deck position.
- **Emergency hand pump:**
allows to use the lift in safety conditions in case of electric installation failure.
- **Gear-motor disengagement:**
for emergency manual manoeuvres.
- **Power supply general disconnecting switch (battery isolator system).**
- **Reflecting strips bounding platform overall dimensions.**
- **Antiskid lift surface.**
- **Intermittent lights of platform overall dimensions.**

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2.4. Identification Plates and Stickers

An identification tag and safety stickers have been applied to the lift, as follows:

WARNING

In the case of loss or excessive wear, identification tags and safety stickers should be rapidly replaced. Please contact an authorized Mobility Networks Ltd centre for details.

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N° MATRICOLA DEL SOLLEVATORE.

3. STABILITY VERIFICATION

Before starting to position and fix the lift, as described in Chapter 4 it is necessary to verify the stability of the vehicle with the lift installed and with full load.

Find following the detailed instruction for stability verification.

WARNING

Please follow the instructions in this manual carefully when carrying out the stability verification of the "RPS 300" wheelchair lift. If the instructions in this manual have not been properly understood or if further information is required, please contact Mobility Networks Ltd immediately. Serious damage to vehicles and human injury could be caused if these instructions are not carefully fulfilled.

WARNING

First of all installer must fill the Use and Maintenance manual at page II and indicates:

CUSTOMER NAME (OWNER OF THE VEHICLE WHERE THE LIFT IS INSTALLED)

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3.1. Side Installation

- Obtain from the vehicle registration documents, the follow information required and write them in Section A of “Log verification for the first commissioning” of the “Use and Maintenance Manual” after have sign the box relevant of lateral installation:

C = track [mm];

T = vehicle mass in running order [daN].

- Get from **Table 1** using **Fig. 3** as reference, data of the lift to install, find following a synthetic description of these data:

G_1 = load of the firm part of the lift, see table [daN];

G_2 = load of the mobile part of the lift, see table [daN];

P = load of the nominal capacity (300 daN);

$M = G_1 + G_2$ [daN];

d_1 = absolute distance between cassette front and of G_1 , application point, see table [mm];

d_2 = absolute distance between cassette front and of G_2 , application point, see table [mm];

d_p = absolut distance between cassette front and P , application point, see table [mm];

- Calculate the B_{min} quote by means of the following formula:


Fig. 3

$$B_{min} = \frac{[(G_1 \times d_1) - (G_2 \times d_2)] - (P \times d_p)}{M + P} + \frac{T \times C}{2.5 \times (M + P)}$$

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B_{min} is the minimum distance between cassette lift front and external edge of the vehicle bumper [mm], that ensures vehicle stability; return the calculated value in the box of side installation of Section A of "Log verification for the first commissioning" of the "Use and Maintenance Manual".

SOLLEVATORE LIFT		LUNGHEZZA BRACCI ARMS LENGTH	G_1	d_1	G_2	d_2	P	d_p
	PORTATA CAPACITY							
	[kg]							
		L [mm]	[daN]	[mm]	[daN]	[mm]	[daN]	[mm]
RPS	300	650					300	

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4. LIFT POSITIONING AND FIXING

In this chapter there are all the necessary instructions to lift "RPS 300" wheelchairs lift on most vehicles on the market. Since it's not possible to cover all possible cases we advise you to contact **Mobility Networks Ltd** for any further information.

Before operating, check that you have the "Application study" document made by **Mobility Networks Ltd** technical department (find following **Fig. 4** is example of application studies for side installation; verify also that it concerns the vehicle you are operating on and the lift, you're going to install.

WARNING

Without this document, you do not have the required complete installation instructions for the correct operation of the lift.

The electro-hydraulic lift warranty becomes void when the assembly has been incorrectly carried out and/or when it does not comply with Mobility Networks Ltd instructions.


Fig. 4

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4.1. Positioning and Verification

4.1.1. Undervehicle chassis installation (lateral)

For this installation a trolley suitable for raising the empty lift with adjustable height deck, is necessary: trolley must be able to bring lift cassette at the distance from the ground listed in the application study (see Fig. 4); then proceed as follow:

- Fix firmly the lift on the trolley deck tacking care not to cover cassette sides fixing plates;
- Open the vehicle doors where lift is going to be installed;
- Adjust trolley deck height in order to bring cassette lower edge at the some quote listed in the application study;
- Put fixing brackets on the lift (supplied in equipment or however approved by **Mobility Networks Ltd**) at the body vehicle;

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Verificare che / Check that:

- Following the instructions of *Chapter 6* of “Use and Manitenance Manual”, open manually the lift and play the manual pump until platform reaches lifting end stroke;

WARNING

Lifting adjustment of the lift is made Mobility Networks Ltd ; check however that the rising quote is in compliance with the application study (take care that internal flap lower must be approximately at $20 \div 30$ from vehicle deck);

WARNING

If quote indicated in application study differ from the real situation of the vehicle on which the installation is being carried out, please, reprogramming lifting angle (see sheet 8.8) or contact Mobility Networks Ltd

- Push the trolley under the vehicle keeping the platform at the outer of the vehicle, enter the trolley internal flap lower edge overcome vehicle deck for 30 mm, then brake the trolley;
- Measure the obtained “B” quote (in compliance with *Chapter 3*); “B” is the distance between cassette front and external edge of vehicle and wheel medium axis (for side installation);
- And check that: $B \geq B_{min}$;
- Using the obtained “B” value calculate the maximun working load of the lift by means of the following formula:

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- Return "B" quote R value calculate in Section A of the "Log verification for the first commissioning" of the "Use and Maintenance Manual", in the section corresponding to the installation you are carrying out (rear or side);
- Check that the fixing brackets are located by vehicle frame or structural elements (or suitable mechanical resistance elements);

WARNING

If the inspection gives negative results, stop the installation procedure. An unstable and unsteady lift installation damages lift correct operation as well as users safety.

- Using fixing brackets as drilling templates, mark the centers of holes to be made on the vehicle structure elements;
- Following the indications of Chapter 6 of "Use and Maintenance Manual", remove the lift from the vehicle and closed it manually.


Fig. 5

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
4.2. Lift Fixing

- Mark the centers holes to be made on the vehicle structure elements; these holes can be:
 - If dead, M8 threaded (for M8 bolts) or M12 threaded (for M12 bolts),
 - If passing, diameter $\varnothing$ 12 mm (for M10 bolts) or diameter $\varnothing$ 14,5 mm (for M12 bolts);
- Replace the lift in the boot space by centering the fixing brackets holes with the ones just made and, by using a suitable length M10 (o M12) bolts, operate the fixing in each plates holes and close them with 52 Nm (or 91 Nm for M12 bolts) clamping torque.

WARNING

If it is not possible to operate the standard installation because it is necessary to have special fixing brackets, follow Mobility Networks Ltd special instructions which have been arranged with the vehicle Installation Centre.

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Example for fixing brackets for FIAT DUCATO 2006

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5. HYDRAULIC SYSTEM

In order to assure to the electro-pump "2" (Fig. 6) an optimally flow of hydraulic oil, it is necessary that the tank "3" (Fig. 6) is always full of oil.

It is possible to check the hydraulic oil level into the tank by unscrewing and removing the cap "7" (Fig. 6) with the platform inside. If it is necessary to fill up oil "**TEXACO RANDO HDZ 15**", check the level oil in the tank is above the "MIN" mark and under the "MAX" mark. At the end of the operation, screw the cap.

⚠WARNING

When the vehicle engine is switched off, it is suggested do not to operate the hydraulic power unit for more than one minute at a time, in order do not to overload batteries.

⚠WARNING

When filling up oil, do not mix oils having different features and gradations.

⚠WARNING

When checking and filling up oil, keep platform on the ground.


Fig. 6

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5.1. Hydraulic Connection

Place in the area in charge of the vehicle the hydraulic unit and then fix it.

To connect the hydraulic system installed into the lift to the external hydraulic power unit, it is enough to connect the flexible hose quick connection coming out from the lift structure to a socket "Z" supplied on the power unit: check that the flexible pipe coupling is tightly hooked, in order to avoid its release when the installation is pressurized (you should hear a click during the insertion, and after the connector should not release if you take out it).

⚠WARNING

It is forbidden to use additional equipments, such as shifting spanners or other objects, because the hydraulic connection between the hydraulic power unit and the lift internal hydraulic installation occurs by means of a swift pressurized connection.

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5.2. Pressure Socket

WARNING

These instructions are for the exclusive use of appropriately trained technical personnel.

DANGER

Before using the pressure regulator make sure that platform is firmly on the ground in order to eliminate the pressure set at zero from hydraulic system pressure.

Inside the hydraulic control unit there is a pressure socket. To measure the pressure in the hydraulic circuit, use an allen wrench 6 to screw off the lid of the pressure regulator placed on the T-piece found on the top of the distribution block of the control unit, in correspondence to the pressure regulator and connect through the male connector 1/4".

WARNING

After disconnecting from the pressure socket, make sure to screw and tighten the lid.

Mobility Networks Ltd decline all responsibility for damage caused by non fulfilment of these instructions and automatically deem the warranty to be void in such cases.

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6. WIRING

6.1. Electrical Connections

Electronic power unit control box is normally located in the vehicle internal part near lift exit door, or inside the vehicle boot sheltered from water and dampness.

WARNING During electrical connections operation, make sure not to have any cables in contact with heat sources and sharp edges in order to avoid damages.

To make electrical connections between lift and hydraulic unit:

Machine wiring exits from the structure right side in rear position and it is contained in a flexible wrinkled special gaiter inside the gaiter there are the following laid down elements:

- Hydraulic power unit wiring ending up with a 6 poles connector (Fig.7) connect it at the connector on hydraulic power unit;
- Driver remote control (keypad) wiring ending with rectangular 35-poles connector (Pos. 1 - Fig.8), connect it on electronic power unit.


Fig. 8

6.2. Electronic control box

- 1) Socket for push-button panel remote control panel;
- 2) "EXIT/LIFTING" and "LOWERING" emergency selector;
- 3) "RETRACTING" emergency selector;
- 4) "RESET" button;
- 5) Controls selector from control panel (emergency) or from push-button panel;
- 6) Working power pack light.

Electronic control box (Fig.8) is normally located in the vehicle internal part near lift exit door, or inside the vehicle boot sheltered from water and dampness.

Selectors "2" and "3" are used only in case of emergency, when happens a keypad failure.

DANGER - Do not leave the lift unattended with the control enable-key inserted in the control panel switch.

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6.3. Vehicle Communication

6.3.1. Lift not stowed

Lift is set in order to send to the vehicle a signal when it is not stowed by means of 2 PINS WHITE FAST-In / FAST-On MALE connector (on lift wiring): use a 2 PINS FAST-In / FAST-On FEMALE connector (for 6.3 male contacts) to connect vehicle to the output signal.

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6.4. Lift Power

- Mount in vehicle batteries compartment the 100A fuse, in the scope of the delivery, connect to one of its poles one 25 mm² red wire and lay this last to reach battery isolator key pole;
- By means of the ø8 lug, in the scope of the delivery, connect the red wire to battery isolator key pole and cover it by means of the rubber cup, in the scope of the delivery;
- Connect the other pole of the 100 A fuse to vehicle batteries positive pole by means of a 25 mm² wire;
- By means of one 25 mm² black wire connect to the ground (-24V) powerpack electric motor negative pole (marked with "-").

 **WARNING:** The connection to the vehicle batteries must be carried out by passing towards the vehicle electric battery-isolator.

 **WARNING:** It is suggested to protect the frame-cable lug with a special paint for electrical protections.

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7. VERIFICATION FOR FIRST COMMISSIONING

The following instructions are necessary for the first commissioning of the **"RPS 300"** wheelchair lift after it has been installed on the vehicle. Note that every record is equipped with a description of the verification check, time estimations, number of personnel required and necessary materials and equipment.

After each verification has been carried out, complete, sign and stamp *Section B* of *Chapter 1* of "Use and Maintenance" manual ("Log verification for the first commissioning").

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WARNING

For the verification check of first commissioning of the **"RPS 300"** wheelchair lift the instructions detailed in this manual must be carefully fulfilled.

If the instructions in this manual have not been properly understood or if further information is required, please contact Mobility Networks Ltd immediately. Serious damage to vehicles and human injury could be caused if these instructions are not fulfilled.

Mobility Networks Ltd declines all responsibility if the verification check for the first commissioning is not followed properly and if the appropriate records are not compiled correctly (in *Section B* of *Chapter 1* of "Use and Maintenance" manual) and will automatically consider the warranty void.

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NOTE

To carry out the static deformation test, follow the instructions below:

- Position the platform at all out position with lifting arms horizontal.
- Measure the height of the platform and its angular alignment in relation to the loading area of the vehicle.
- Apply a load of 375 kg onto the platform and then remove it.
- Repeat the height and angular measurements of the platform and check that there are no permanent deformations either to the lift or the vehicle fixings which could be a detriment to the correct functioning of the lift.

To carry out the static derivation test, follow the instructions below:

- Position the platform at the loading area level of the vehicle.
- Apply a load of 375 kg onto the platform.
- Measure the height of the platform and its angular alignment in relation to the loading area of the vehicle.
- Repeat the measurements after a trial period of 15 minutes. Check that the vertical derivation of the platform between the two measurements is not greater than 15 mm and that the angular derivation is not greater than 2°.

To carry out the dynamic test, follow the instructions below:

- Position the platform at ground level.
- Apply a load of 300 kg onto the platform.
- Check the complete functioning of the lift's up and down movements when fully loaded.

See Annex II of EN 1756 rule

To carry out the safety check of over load:

- Position the platform at ground level.
- Apply a load of 375 kg onto the platform.

The check is positive if the platform don't leave ground level.

If check is negative contact **Mobility Networks Ltd**

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8. STANDARD INSPECTIONS

Standard inspections should be carried out 6 month by **Mobility Networks Ltd** Authorized Service Centre. They allow for greater operational safety and guarantee a longer lifetime for the **RPS 300** wheelchair lift.

Here are the instructions for standard lift inspections to be followed by service centre. In particular, in every inspection form you can find a description of the verification carried out, and of the possible intervention in case of a negative check test result, time estimations, number of personnel required and necessary materials and equipment.

The record for the standard inspection, present "Use and Maintenance" manual, must be filled, stamped and signed after every test check.

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WARNING

In the standard inspection phase, the servicer is also required to carry out any maintenance interventions even if they are the responsibility of the user of the lift, and as described in *Chapter.7.* of "Use and Maintenance" manual.

WARNING

After **2000 HOURS** of functioning and always after **10 YEARS**, the lift should undergo a full and complete service at a **Mobility Networks Ltd** authorized centre, who will subsequently issue a new user manual with new inspection procedures prior to resuming service.

IF THIS IS NOT THE CASE THE LIFT MUST NOT BE OPERATED.

WARNING

For the standard inspections on behalf of the servicer of the **RPS 300** wheelchair lift the instructions detailed in this manual must be carefully adhered to.

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Mobility Networks Ltd decline all responsibility for damage caused due to incorrect maintenance or lack of maintenance of the lift and automatically deem the warranty to be void in such cases.

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Park the vehicle on a level and horizontal (as much as possible) surface , by taking care that any objects are located on it and that this surface is clean. Here, you can carry out the lowering operation.

“Extension limit position” adjustment

Always check that the platform is completely exit from its housing.

Otherwise, it is necessary to adjust the lock distance of the micro “all outside” to do this operation, fulfill the following indications:

- Unscrew the lock screw,
- Move the lock, of the quantity required, by following the special slot direction,
- Close the screw.

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Cylinders running adjustment (loading level height)

Always check that the platform lifts as far reaching the loading level.

However, if the platform does not reach the loading level, it is necessary to adjust the locks of both cylinders reducing the movement. In order to do that, fulfill the following indications:

- Unscrew the "G" grains with an Allen spanner.
- By using an open spanner, rotate the cylinder head.
- Once having adjusted the running, tighten the grains again.

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PLATFORM REACHING TO THE HORIZONTAL PLANE

When lifting upper arms are horizontal, check with a blister that the platform is horizontal (vedi figura a lato).

If the platform is not on the horizontal plane, it is necessary to record the lifting lower arms length, by loosening the stop nuts and rotating the stud bolt.

Once having carried out the operation on both lower arms, start stopping the stop nuts again.

WARNING

By screwing the stud bolt, the arm lengthens; otherwise, by unscrewing it, the arm shortens.

WARNING

The length of the two lower arms must be identical: take care during control.

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10. ALLEGATI / TECHNICAL ANNEXES

9. ADDITIONAL INSPECTIONS

The authorized **Mobility Networks Ltd**, service centre which has been contacted by a customer encountering anomalies or malfunctioning of the wheelchair lift “**RPS 300**”, must carry out the required checks and if necessary operates in order to solve the associated problems. At the end of the additional inspection service centre must fill, stamp and sign the relevant form in lift “Use and Maintenance” manual.

WARNING

In the additional inspection phase, the service centre is also required to carry out any maintenance interventions even if they are the responsibility of the user of the lift, and as described in *Chap. 7* of Use and Maintenance manual.

WARNING

For the additional inspections on behalf of the servicer of the “**RPS 300**” wheelchair lift the instructions detailed in this manual must be carefully adhered to.




If the instructions in this manual have not been properly understood or if further information is required, please contact **Mobility Networks Ltd** immediately. Serious damage to vehicles and human injury could be caused if these instructions are not adhered to.

Mobility Networks Ltd decline all responsibility for damage caused due to incorrect maintenance or lack of maintenance of the lift and automatically deem the warranty to be void in such cases

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
COLLEGAMENTI INTERNO CARRELLO / WIRING CONNECTIONS ON TROLLEY

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

CONNETTORE 35 POLI SEALED

WIRING POSITION	
1	LIGHT PLATFORM
2	ELECTRIC VALVE
3	TELERUPTOR
8	GEAR MOTOR
9	GEAR MOTOR
10	GEAR MOTOR
11	GEAR MOTOR
12	ELECTRIC ACTUATOR
17	CABIN LIGHT
14	BATTERY
18	STOWING MICRO
19	PRESSURE SWITCH
20	MICRO TD/TF
21	PLATFORM STOWING
22	DEPLOY PLATFORM
15	LIGHT MICRO
16	LIGHT MICRO
24	BATTERY +
25	BATTERY +
26	BATTERY -
27	BATTERY -
28	BATTERY -
31	SENSITIVE EDGE
33	BATTERY +
35	ELECTRIC ACTUATOR

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- 1 FILTER 90 MICRON
- 2 ELECTRIC PUMP 12V dc 500W
- 3 GEAR PUMP 0,48 cc/g
- 4 MAX PRESSURE VALVE 14MPa
- 5 ELECTRIC LOWERING VALVE
- 6 EMERGENCY MANUAL PUMP
- 7 CHECK VALVE
- 8 LOWERING MANUAL VALVE
- 9 ADJUSTABLE VALVE
- 10 LIFT CYLINDER

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


***REFER TO MOUNTAIN DIAGRAM FOR THE BRACKETS POSITION.
FIRST OPERATION : MOUNTAIN THE BRACKETS ON LIFT.***

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CHECK THAT THE PLACEMENT OF THE LIFT CORRESPONDS OPENING DOORS OF THE VEHICLE, INDICATIVELY THE INDICATE MEASURES ARE BY CHECKING DURING THE MOUNTING OF THE LIFT

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**USING A CRANE TO LIFT THE VEHICLE
AND DISMOUNTING THE SILENCER TUBE**

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
FOR MOUNTING THE LIFT ON THE VEHICLE, FIX THE REAR RIGHT BRACKET BECAUSE IT IS SLOT (B) IS REFER TO ASSEMBLY

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



INSERT THE DISTANCE IN THE HOLES ($\varnothing$ 18) AND FIXING THE BRACKET ON THE VEHICLE FRAME HOW TO APPLICATION STUDY.

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



Mobility Networks Holdings Ltd

12 Estuary View Business Park, Whitstable, Kent, CT5 3SE, UK

Tel: +44 (0)1227 505022 - Fax: +44 (0)1227 505023 - www.mobilitynetworksgroup.com

Mobility Networks is a trading name of Mobility Networks (Holdings) Ltd.

Registered in England & Wales number 08295602